

\$750

PhaSeal® Starter Kit

Injection Safety Products *Injector Luer Lock*



Encapsulated cannula that attaches to a disposable syringe with a Luer Lock connection.

PRODUCT CODE

QUANTITY

N31

25

Drug Safety Products



Protector 14 For drug vials with a 14 mm neck.
Use with syringes up to 20 ml.

P14

6

Protector 50 For drug vials with a 21 mm neck.
Use with syringes up to 50 ml.

P50

6

Protector 21 For drug vials with a 20 mm neck.
Use with syringes up to 20 ml.

P21

6

Protector 28 For drug vials with a 28 mm neck.
Use with syringes up to 60 ml.

P28

2

Administration Safety Products *Connector Luer Lock*



Attaches to the patient's IV line, providing a closed connection to injector.

C40

10

Y-Site



Y-site Connector provides a dry connection port for an IV push without interfering with other ports that are in use

C80

10

Infusion Adaptor



Connection between infusion bag and fusion set, with built-in Connector.

C100

10

Application Products *Assembly Fixture*



Mechanical device that facilitates the attachment of Protector to drug vial.

M12

1

Pricing Also Includes:

Clinical Training, Educational and Regulatory Information, and Continued Clinical Support



Dr. Brett Cordes, DVM
Medical Director – Veterinary Pharmacy

866.602.4707 Toll-Free Phone

888.602.2242 Toll-Free Fax

bcordes@theapothecaryshop.com

PhaSeal®

theapothecaryshop.com